

10 LITRES

Product Code : MAOXY10 / (AE-8-S)



- Flow rate 1-10 LPM, with an increment of 0.5 LPM
- Can be used for pediatric also
- Provides a continuous flow of concentrated oxygen (93% ± 3%) from room air through one oxygen outlet at maximum flow of 10 LPM
- PSA (Pressure Swing Adsorption) Technology
- Inbuilt Nebulizer and Timer (10 MIN to 5 HOURS)
- Alarms for Low Oxygen purity, Low pressure, Power Outage, High/Low voltage alarm and defective information is displayed on the screen
- HEPA filter and Bacteria filter. Five-stage filters keep away most of impurity, bacteria and PM (particulate matter) 2.5 microns.
- Run with Stabilizer (2KVA, 6AMPS)

Specifications :

Flow rate	1-10L/MIN
Concentration	93% (± 3%)
Power Consumption	610VA
Sound Level	<60db
Outlet Pressure	100 Kpa ±10% or 14.5 psi±10%
Power failure Alarm	Yes
Over-heat Alarm	Yes
Low Purity Alarm	Yes
Voltage Alarm	Yes
Low Pressure/High Pressure alarm	Optional
Maintenance prompt	Yes
Power On self test and LVD	Yes
Designed for 24 hour working	Yes
Net Weight	24 kg
Gross Weight	26 kg
Warranty	2 years / 12000 hours



Years Warranty*



Prompt and Efficient
Service



Dual Filtration



EN ISO
13485:2016