

Karma



Aqua 2 CMD-MS Blue

New range in commode chair with best quality material and 1st in segment commode with **120KG** weight bearing capacity

CE TÜV Rheinland ISO 13485

AQUA
COMMUNE SERIES

1 YEAR
WARRANTY

terms & conditions apply

The images shown are for illustration purposes only and may not be an exact representation of the product

Aqua 2 CMD-MS Blue

Karma Aqua series commode chair made from MS with partial elliptical frame design. Aqua series commode follow the best industrial standard powder coating which provide best protection and environment friendly.



Perfect Seating Comfort
unleash western seat
comfort in style



Injection Molding backrest
provide extra durability



Revolutionary Ergo Bucket
design redefining ease of
use



**Foldable frame ease for
storage**



Height Adjustable leg
provide convenience as
per user height

TECHNICAL SPECIFICATION

Frame Style	Frame Material	Total Width	Seat Width	Wheel to Wheel Width (Close)	Rear Wheel Size	Front Wheel Size	Seat Height Adjustment (25mm adjustment)	Seat Depth	Back height	Total Length	Max user Weight Capacity	Net Weight	Gross Weight
Standard	MS	532 mm	460 mm	340 mm	3"	3"	773-873 mm	426 mm	323 mm	535 mm	120 kg	7.87 Kg	8.90 Kg

measurement can be deviate +5mm #Karma reserve all rights to change product information without any prior notice



Marketed by
karma
Healthcare

KB 22, Bhakta tower, 9th Floor
Sec 3, Salt Lake, Kolkata 700098

The images shown are for illustration purposes only and may not be an exact representation of the product