



**-Life®**



**24 Hours**  
Continuous working



**High Purity**  
Oxygen 93%



**Wise Alarm**  
Function



**Five-channel**  
Oxygen Filter

## OC-102 5 Litre Oxygen Concentrator

### Features:

- Running Time Calculation
- Protection for over-heat & over load
- Power failure alarm
- Alarm for high/low pressure
- Purity alarm when the purity below 82%
- Medical Nebulization
- Low flow alarm
- High temperature alarm

### Specifications:

Oxygen Flow	1-5 L/minutes @ 93% ±3%
Output Pressure	6 - 10 PSI
Purity	±3%
Comparative humidity scale	≤95%
Power Consumption Noise	≤45 dB(A)
Oxygen Concentration	93% ±3%

### How to Use:



Fill the humidifier with water upto the max level.



Easy installation with hands.



Fix the breathing tube.



Switch on the button to start the oxygen flow.



Use LPM knob to maintain the oxygen flow.

### Contact us:



Mobile No./ WhatsApp No. : +91 9310113243



Web : [www.klifecare.com](http://www.klifecare.com)



Email : [info@klifecare.com](mailto:info@klifecare.com)

**Kannu Impex (India) Pvt. Ltd.**

412, Suncity Trade Tower, Sector-21, Gurugram-122016