



Ryder 1

LESS IS MORE



Powder coating frame



Comfortable position of armrest



Attendant brake included



Soft and breathable upholstery & Safety belt



Fixed Footrest



16" solid rear wheel & 8" solid castor

Seat Width	Seat Depth	Overall Width	Fold Width	Armrest Height	Overall Length	Seat Height	Backrest Height	Overall Height	Weight	Total Weight
18"	16"	24"	9"	10"	33"	18"	17"	34"	13Kg	100 Kg

Karma Healthcare Ltd. KB 22, Bhakta Tower, Sec 3, Saltlake Kolkata 700 098, Customer care ☎:033 2335 7062 ✉:service@karmahealthcare.com