




RYDER 14

aluminum wheelchair



Traveling



Home use



Hospital use

FEATURES



FLIP BACK ARMREST

Armrest can flip back make the user transfer easy

FOLDABLE BACKREST

Wheelchair backrest is foldable which make easy to store



FLIP UP FOOTREST

Footrest can flip up which help to move in tight space

ATTENDANT BRAKE

The model is equipped with attendant brake for safety



TECHNICAL SPECIFICATION

Frame Style	Frame Material	Wheel to Wheel Width(open)	Seat Width	Wheel to Wheel Width(Close)	Rear Wheel Size	Front Wheel Size	Seat to Floor Height In	Seat Depth	Back height	Total Height	Total Length	Max user Weight Capacity	Net Weight	Gross Weight
Standard	Aluminum	24 inch	18 inch	12 inch	22"	6"	19.5 inch	16 inch	16 inch	35 inch	41 inch	100 kg	12.50 kg	14.67 Kg