



Ryder 5

Manual Wheelchair



Detachable Footrest
Convenient for movement in
tight place



Foldable Backrest
ensure travel friendly
mobility



Flip back Armrest
Help for easy user
movement



Improve quality
Upholstery for comfort



Attendant Brake for
safety



Compact and Foldable
can take any where

Seat Width	Seat Depth	Overall Width	Fold Width	Armrest Height	Overall Length	Seat Height	Backrest Height	Overall Height	Weight	Total Weight
18"	16"	25"	10"	10"	36"	18"	16"	32"	16Kg	100 Kg

Marketed by: Karma Healthcare Ltd. KB 22, Bhakta Tower, Sec 3, Saltlake Kolkata 700 098