



## Ryder MWC 6 F12



New range of manual wheelchair in MS range. Light weight & comfortable. Best for indoor use

# Ryder MWC 6 F12

We are pleased to introduce our latest mild steel attendant-driven wheelchair. This wheelchair is designed to be compact and suitable for indoor use, while also being perfect for travel due to its portability. It is equipped with the best features for your convenience.



**Flip up footrest** makes the wheelchair more compact



**Fold able frame** and compact for easy transportation



Extra safety for user comes with **attendant packing brake**



**Backrest** of the wheelchair can also **fold able** for easy storage



**Flip up armrest** makes easy for user shifting from one place to another

## TECHNICAL SPECIFICATION

Frame Style	Frame Material	Wheel to Wheel Width(open)	Seat Width	Wheel to Wheel Width(Close)	Rear Wheel Size	Front Wheel Size	Seat to Floor Height In	Seat Depth	Back height	Total Height	Total Length	Max user Weight Capacity	Net Weight	Gross Weight
Partial Elliptical	MS	584 mm	457 mm	241 mm	12"	6"	483 mm	394 mm	381 mm	851 mm	858 mm	100 kg	11.57 Kg	13.80 Kg

# measurement can be deviate +-5mm #Karma reserve all rights to change product information without any prior notice



Marketed by  
**karma**  
 Healthcare

KB 22, Bhakta tower, 9th Floor  
 Sec 3, Salt Lake, Kolkata 700098